EVALUATION OF DYNAMIC WAVEFRONT ABERRATION MEASUREMENT

DANIEL NEAL, PD, XIFENG XIAO, PHD, TOM DUNLAY, MATT HAUGO

WAVEFRONT DYNAMICS INC...

Purpose

The eve is in constant motion, but most objective measurements are taken as a single measurement. By recording a sequence of measurements, a better evaluation of the eye condition can be obtained. This is particularly important for correcting aberrated eyes with custom contact lenses, where translation and rotation of the correction directly affects the visual quality. In addition, accommodation may play a significant role in the calculation of refraction from a wavefront measurement.

There is often a difference between a subjective manifest refraction and objective measurement.1,2 Autorefractors and aberrometers both suffer from "instrument accommodation" where the eye accommodates during the measurement process. There are several factors that are important:


Fixation target

- Image content
- Astigmatism correction
- Fogging
- Eye condition
 - Tear film
 - . History
 - Refractive state

These factors make it complicated to determine basic refraction

Most autorefractors and aberrometers measure the eve with a single snapshot. By recording a dynamic sequence³, our hypothesis is that we can improve on the objective refraction accuracy.

TYPICAL DYNAMICS



Dynamic refraction parameters

Methods


DYNAMIC ABERROMETER

A dynamic wavefront measurement instrument was developed using a high-resolution Shack-Hartmann wavefront sensor. This was combined with a large field of view imaging camera and corneal topographer to provide a platform for dynamic measurement. The system used two synchronized cameras to measure the eye at up to 54 fps. Instrument accommodation was studied by comparing maximum measured sphere to average sphere in a 10 second sequence, for five subjects, ages 22-42. Subjects were measured five times each with three different operators.



Specifications


- A multi-function system capable of
- measuring Measurable range: sphere and Dynamic wavefront aberrometry cylinder measurements in 0.01 D Dynamic iris image
- Contact lens fit increments. Spherical equivalent range (7 mm pupil) -22 to +16 D Refraction Cylinder range (7 mm pupil) 14 D Dynamic corneal topography
- Axis in 1° increments Dynamic keratometry Pupil dynamics Maximum wavefront diameter 81
- Accommodation range 24 mm horizontal field of view Zernike terms displayed through
- sixth orde Measurement spatial resolution 116 µm (approximately 2870
- measurement points for a 7 mm pupil
- Eye image field of view 24 mm
- Subjective refinement with integrated evechart
- Live wavefront and refractive displays
- Automated patient alignment stage



Results and Discussions A

DYNAMIC ABERROMETER

For all subjects, the measured sphere varied during the 10 second measurement sequence. Data was acquired at 10 or 54 fps




rs (Seq. Sph and Cvl) as a function of time for Refraction pa a 23 y.o. female subject



FINDING OPTIMUM REFRACTION

- · Merit function parameters
 - · Seq, pupil size, image quality
 - Find maximum MF in sequence
 - · Adjust weights to match clinical data

$MF = W_{SEQ}SEQ + W_PP_D + W_{IQ}IQ$










Refractive parameters with merit function calculation. W_{seq} =3, W_{ρ} =1, W_{lq} =1. The maximum MF is found with a simple search.

DATA ANALYSIS









The maximum (most hyperopic) measurement was on average 0.25D more than the average sphere. In the young group (ages 20-25) the average difference was 0.28 +/-0.1 D while in the old group (age 38-42) the difference was 0.21+/-0.07 D. The max sphere values were more myopic by -0.34+/-0.35D compared to a commercial aberrometer

Conclusions

totio bentue totion bentue bentue bentue bentue bentue

Outp Institute Coup Institute Coup Institute Institute Institute Institute

The MF selection matched closer to iDesign refraction values than the sequence average. Overall, compared to iDesign, the NextWave measurement was slightly myopic in the younger group but hyperopic in the older. This may be due to the very simple fixation target in this early prototype which is planned to be improved in future versions

The merit function approach gave more consistent results with fewer outliers.

m Vis Sci 75(8) 617-22 1998

