

Prospects for practical higher order aberration correction

DANIEL NEAL, PHD\*, CHRISTINE SINDT, OD\*,XIFENG XIAO, PHD\*  
\*WAVEFRONT DYNAMICS INC., \*UNIVERSITY OF IOWA

Purpose

About 14% of the population has elevated higher order aberrations (HOA) which result in some loss of visual quality [1]. About 1% suffers from keratoconus, pellucid marginal degeneration, corneal transplant, post-op ectasia or other condition that significantly degrades their vision [2]. Typically, these individuals are not eligible for LASIK and existing modalities (GP lenses, etc.) may not sufficiently restore normal vision. Thus, a measurement driven custom wavefront correction is needed [3]. Many of these patients are fitted with scleral contact lenses. However, these lenses often do not restore vision to normal levels due to incomplete index matching, corneal back surface and lens.

Previous studies have demonstrated the correction of wavefront error with scleral lenses and have shown the effectiveness of this treatment [4-6].


The purpose of this work is to evaluate practical ways to ameliorate barriers to adoption of Higher Order Aberration (HOA) correcting contact lenses (CL).

HOA contact lenses are not new, and several researchers have demonstrated effective HOA correction is possible. Clinicians have expressed reluctance in adoption, citing concerns such as effectiveness, feasibility, equipment, patient variables and cost.

Perceptions of HOA contact lens practicality – social media posts

- Questions extracted from various social media posts
- Are aberrations stable over time?
  - Can we measure the eye with large aberrations?
  - How do we get consistent, accurate measurement?
  - Can we measure through a scleral lens?
  - Are contact lenses stable in position and rotation?
  - Can we make contact lenses consistently and accurately with the desired correction?
  - What patients should we select for HOA correction?
  - How does pupil size and shift affect aberration correction?
  - Does the wavefront correction add steps to the process?
  - Does it work?

Scleral lenses alone do not always improve the net aberration



Methods

Through social media posts, on-line chat rooms, and focus groups of CL practitioners, perceived HOA process issues were identified. Concerns were categorized as: measurement accuracy, workflow, patient factors, and outcomes. We evaluated and optimized factors affecting these concerns. A clinical study using these methods looked at chair time and visual outcomes.

40 consecutive eyes from 20 patients, 54% male, 46% female ages 23-76, were fitted with rotationally stable impression-based scleral contact lenses [8]. A dynamic, high resolution aberrometer measurement was made at each patient visit. The design integrated low order aberrations (LOA) and HOA into the optical zone.

**DYNAMIC ABERROMETER**

A high resolution, dynamic wavefront aberrometer system was developed to provide accurate measurements of highly aberrated eyes [7]. In a feasibility clinical study, this instrument was used to guide the treatment of patients with aberrated corneas. The instrument was used to evaluate the aberrations, provide parameters for the design of the scleral lens, evaluate the fit of the lens, and then a second lens was designed to correct the wavefront aberrations based on measurements through the first lens. For those patients that already were fit with a scleral lens, their existing lens was used as the predicate.

Results and Discussion

**Measurement accuracy:** Measurement accuracy was facilitated by using a high resolution lenslet system, sequential data capture over time through blinks to account for tear film variability, accommodative control in the system to prevent instrument myopia, a high-resolution camera to identify ocular features, and topography for patch centration.

**Workflow** With objective fitting techniques, only 3 office visits were needed: 1) base lens fitting visits 2) dispensing of LOA and HOA measurement visit 3) dispensing of HOA lens (reduced from 5 to 3 visits). LOA and HOA were corrected simultaneously reducing lens remakes from 4 to 2.

**Patient Factors** HOAs are stable over multiple timecales. Dilating the pupil increases the treatment area is bigger than the natural pupil.

**Outcomes** The HORMS reduced from an average of 0.9 to 0.35  $\mu$ m in one iteration, which was a 62%+13% improvement. All eyes had reduced HOA, no eye lost any lines of acuity.

- Dynamic wavefront aberrometry, iris image, pupil, corneal topography, keratometry and refraction
- Contact lens fitting
- Subjective refinement with integrated eyechart
- Live wavefront and refractive displays




Figure 3. Multi-function aberrometer/topographer

Results and Discussions A




Figure 4. WFG lens reduction in aberrations. a) bare eye, b) traditional scleral, c) WF corrected




Figure 5. Long term HOA lens stability on the eye




Figure 6. Long term HOA lens stability on eye a) initial, b) after 6 months of wear

Results and Discussions B

- Workflow**
- Visit 1
- Evaluate LOA scleral lens fit
  - Screen for HOA correction
  - Dilate eye
    - Measure aberration through LOA lens
  - Order HOA lens
  - Dispense LOA lens
- Visit 2
- Evaluate HOA lens performance
  - Dispense
- Clinical study**

The severity ratio for each eye is calculated by dividing the Sph/Cyl corrected HORMS value by the age/pupil normal higher order RMS for that eye. The severity ratio can be used as a guide to predetermine which patients will benefit significantly from HOA correction. In this study 50% of the patients were at least 2X the normal range.




Figure 7. Severity Ratio for the study patients.




Figure 8. HOA RMS for bare eye, scleral lenses and WFG lenses



Figure 9. HOA RMS results for bare eye, scleral lenses and WFG lenses

Results and Discussions C

These scleral lenses are extremely stable on the eye, so the correction works very well in nearly all cases. WFG correction brought the subject to normal aberration levels in most cases.

The WFG refraction also improved the refraction compared to the standard scleral lens since the low orders are included in the WFG correction.

On average scleral lenses reduced the bare eye higher order aberrations (6 mm pupil) from 1.70 $\pm$ 0.87  $\mu$ m RMS to 0.94 $\pm$ 0.33  $\mu$ m RMS (reduction of 52% $\pm$ 16%). The wavefront correction reduced the aberrations to 0.34 $\pm$ 0.19  $\mu$ m RMS (83% $\pm$ 16% reduction). This resulted in a 62% $\pm$ 13% improvement over the basic scleral lens. All patients preferred the wavefront guided lens.

Lens	Percent Improvement
Sph/Cyl scleral lens compared to bare eye	52 $\pm$ 16%
WFG lens compared to bare eye	83 $\pm$ 16%
WFG lens compared to sph/cyl scleral	62 $\pm$ 13%

Table 2 – The RMS comparison of bare eyes, Sph/Cyl lenses and WFG lenses

Conclusions

While there are barriers to HOA correction with contact lenses, these barriers can be identified and resolved. For the current 40 eye study all patients experienced a reduction in HORMS with a single iteration. This provides practical HOA treatment for aberrated eyes for keratoconus and other conditions.

References

1. D. Godeffroy et al., "Age-specific incidence and prevalence of keratoconus: a nationwide registration study," *Am. J. Ophthalmol.* 195, 169-172 (2019).
2. Zadnik et al., "Baseline findings in the collaborative longitudinal evaluation of keratoconus (CLEK) study," *Invest. Ophthalmol. Vis. Sci.* 39(23), 2537-2546 (1998).
3. Hastings, G.D., Asanpette, M.S., Nguyen, L.C., Ho-Pham, M.J., Hennrich, R.T., Marack, J.D., Comparison of wavefront-guided and best conventional scleral lenses after habitation in eyes with corneal ectasia, *Optometry and Visual Science*, (2018) 95(4):238-247, PMC6450417.
4. Marack, J.D., Parker, K.E., No, Y., et al., On-eye Performance of Custom Wavefront-guided Soft Contact Lenses in a Habitual Soft Lens-wearing Keratoconic Patient, *J. Refract. Surg.* 2007;23:360-4.
5. Soltan, R., Acem, J., Corvino, L., et al., Vision improvement by Correcting Higher-order Aberrations with Customized Soft Contact Lenses in Keratoconic Eyes, *Opt. Lett.* 2007;32:1000-2.
6. Neal et al., "Dynamic aberrometer/topographer designed for clinical measurement and treatment of highly aberrated eyes," *Optical Engineering* 32(8):8, December 2022 • Vol. 61(12).
7. Agapito et al., "Three-dimensional relationship between high-order root-mean-square wavefront error, pupil diameter, and aging," *J. Opt. Soc. Am. A* Vol. 24, No. 3/March 2007.
8. Sheeman JM, Hoffman JM, Zimmerman MD, Ling J, Greiner MA. Indications for Wave Visual Outcomes, and Complications of Custom Imprint 3D Scanned Scleral Contact Lens Use. *Cornea*. 2021 May 1;40(5):596-602. doi: 10.1097/ICO.0000000000002588. PMID: 33252387.
9. Cheng, Kimbrough, Kolbaum, Tobias and Bradley, "Test-retest reliability of clinical Shack-Hartmann measurements," *IOVS*, Jan 2004, Vol. 45, 331-340.