



Scleral Lens with Novel Wavefront Guided Extended Depth of Focus Optics for Visual Improvement in a Presbyopic Patient with Keratoconus

John D. Gelles OD1,2, Nicolas Brown BS 3, Andres Guevara PhD 3, Sydney Krisa BS1, Daniel Chung BS1, Steven A. Greenstein MD1,2, Peter S. Hersh MD1,2

1) Cornea and Laser Eye Institute-Hersh Vision Group, CLEI Center for Keratoconus, Teaneck, NJ

2) Institute of Ophthalmology and Visual Science, Department of Ophthalmology, Rutgers New Jersey Medical School, Newark, NJ

3) OVITZ Corporation - Rochester, NY, USA


RUTGERS
New Jersey Medical School

Introduction


- This case reports on the use of wavefront-guided (wfg) extended depth of focus (EDOF) optics on a scleral lens (SL) for visual improvement in a presbyopic patient with keratoconus (KC).

Case Description

- A 49-year-old female with keratoconus (KC) was fit with a traditional optics SL (tSL) for distance correction.
- A wfgSL was created using a comprehensive system (Ovitz, xWave, Rochester NY) that included a dot matrix on the SL and a wavefront aberrometer with iris and dot registration with direct data transfer.
- Once finalized, a novel wfgedofSL was manufactured. Distance visual acuity (DCLVA), near visual acuity (NCLVA), and total higher-order root mean square (HORMS) with pupil diameter matching were measured with each lens.
- Data was collected after 2 weeks of lens wear and a minimum of 3 hours wear prior to examination.



Zernike Coefficient Comparison



Results

- DCLVA for OD, OS, OU were 20/20, 20/20, 20/20 for tSL, 20/15, 20/15, 20/15+ for wfgSL, and 20/15, 20/20, 20/15+ for wfgedofSL.
- NCLVA was 20/30, 20/30, 20/25 for tSL, 20/30, 20/25-, 20/25 for wfgSL, and 20/20, 20/20, 20/20 for wfgedofSL.
- HORMS for OD and OS was 1.34um and 0.61um for tSL, 0.49um 0.23um for wfgSL, and 0.63um and 0.62um for wfgedofSL. OD values are reported at a pupil diameter of 7.4 mm; OS at 6.3 mm.
- Comparison of tSL to wfgSL HORMS was reduced by 63% and 62%, resulting in a DCLVA improvement of 1 line OD, OS and OU and NCLVA improved of 0.5 lines OS and no improvement OS and OU.
- Comparison of tSL to wfgedofSL, HORMS was reduced by 52% and 0%, resulting in a DCLVA improvement of 1 line OD and OU and 0.5 OS and NCLVA improvement of 2 lines OD and OS and 1 line OU. The increase in HORMS from the wfgSL to the wfgedofSL is expected due to the nature of the EDOF correction.

Conclusions

- When compared to tSL, the wfgedofSL had reduced HORMS and improved visual performance at all distances.
- Future larger prospective studies are required to corroborate this data.

Disclosures

Gelles, JD has received research support, devices, or honoraria from Ovitz. Guevara, A and Brown, N are employed by Ovitz. All other authors have no relevant financial disclosures.



The Cornea and Laser Eye Institute
CLEI Center for Keratoconus